

Reminder

Adresses : Users & Riggers

Equipment affected : Pilotchute « Quick 3-1 » P/N. 2005030103.

REMINDER OF THE PACKAGING :

Depending on the Use and Maintenance Manual Atom Legend

P/N P9009113160, Edition 02 of march 2003 and Edition 02.1 of september 2003.

- Compress the pilotchute by pulling up the net towards the middle coils, then the fabric in the upper coils as indicated by the arrows on the drawing below.
- During the compression sequence, the pilotchute must be firmly held, without any coil sliding.
- For further information, please contact Parachutes de France After Sales Service.
Contact : Mr Philippe SCHOLTES Tel : +33 (0)1 34 32 77 57 // e-mail : pscholtes@zodiac.com

